



Liquid Oxygen System with Conserving Device

Liquid oxygen systems provide excellent mobility for patients who do not want oxygen therapy to compromise their active lifestyle in any way.

The base unit provides patients with enough oxygen for 9 days' continuous flow* and is designed to make receiving oxygen therapy indoors extremely easy. Not only that, but the detachable lightweight portable unit makes it possible for patients to carry on with a normal life beyond the home too.

When used with a conserving device the lifetime of the portable unit can be further increased by between 3 to 5 times, depending upon the flow rate.



Total flexibility...

- **Mobility** – the standard portable unit weighs only 3.6 kg and provides up to 8 hours of ambulatory oxygen whilst away from home. The lightweight system is extremely compact, weighs just 2.5 kg and can provide up to 4 hours continuous oxygen supply*
- **Easy to monitor** – both the base and portable units contain contents gauges to show just how much oxygen is available at any time.
- **Simple to use** – quick and easy to fill; it takes about 40 seconds
- **Efficiency** – can be used with a conserving device to increase the lifetime of the portable unit
- **Flexibility** – a wide range of flow rates are available between 0.25 to 15 l/min
- **Ready to go** – bag, mask cannula and other accessories provided

* based on 2 l/min

We don't take breathing for granted.

Easy to use...

With liquid oxygen systems, patients can return to a more active lifestyle without the encumbrance of bulky, heavy equipment and the bother of frequent refilling.

The liquid portable units offer a range of flow rates between 0.25-15 l/min. It is therefore ideal for patients who require higher or lower flows other than 2 or 4 l/min.

On the initial visit, a Vitalair patient service representative will deliver and install the equipment and provide the patient with comprehensive product training. They will then be responsible for replenishing the patient's liquid oxygen on a regular basis.

The Vitalair team are always on hand to answer your queries and can be contacted 24 hours per day on freephone number **0800 136 603**.



Specifications:

The Companion

	<i>Dewar</i>	<i>Standard</i>	<i>Lightweight Portable Unit</i>
Length	1020 mm	350 mm	270 mm
Diameter	356 mm	155 mm	150 mm
Weight (full)	61.2 kg	3.7 kg	2.5 kg
Gas capacity	25,650 litres	1.2 litres	0.5 litres
Flow rate	0.1-15 l/min	0.25-15 l/min	0.25-6 l/min
Duration of supply (without a conserving device at 2 l/min)	At least 220 hrs	8 hours	4 hours

Actual contents, dimensions and weights may vary from the figures indicated.

**BOC Medical
Vitalair Team**
Customer Service Centre
Priestley Road
Worsley
Manchester
M28 2UT

Tel: 0800 136 603
Fax: 0800 169 9989

BOC Medical is a trading name used by operating companies within The BOC Group, the parent company of which is the BOC Group PLC. The stripe symbol and the words BOC and Vitalair are registered BOC Group trademarks. Copyright The BOC Group PLC 2003.